

Tutorial for EasyPubPlot

PCA/PLS-DA Scores Plot

Step 1: Prepare the Input Data

1. Prepare the input data file with columns of PCs and rows of sample names.

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O
1		PC1	PC2	PC3	PC4	PC5	PC6	PC7	PC8	PC9	PC10	PC11	PC12	PC13	
2	C126	-2.5672	5.5404	1.7257	-1.831	0.67347	-3.3714	2.7031	-0.33016	1.853	-0.08743	2.7316	0.32416	-2.57E-15	
3	C132	-1.4291	5.6798	1.6342	1.3839	1.39	4.4081	-3.3997	-1.8444	1.3685	-0.1636	0.045095	-0.19257	2.30E-15	
4	C152	9.8425	3.2665	-4.0284	-1.2067	2.4641	-0.18291	0.39152	-0.056	-2.082	1.9663	-0.082	-0.53847	-1.31E-14	
5	C153	0.61538	2.8167	-3.7068	-0.54344	-6.9167	0.017108	-0.70333	1.3799	0.7068	-0.74869	-0.65962	0.03615	1.39E-15	
6	C127	-5.8597	1.3255	-0.70983	4.5108	0.89851	0.80902	3.6458	1.1911	-0.75092	0.58059	-1.9261	-1.7402	1.39E-15	
7	C135	-1.1653	1.0327	3.7007	-1.0013	0.73813	-1.2317	-1.0538	2.3779	-1.0516	0.99118	-2.3749	3.4032	1.02E-15	
8	C138	-1.9439	-0.97578	0.8147	1.8415	0.34581	-4.4083	-3.8613	-0.05883	-2.1069	-1.1884	0.60605	-2.0239	-4.66E-16	
9	C160	1.299	-3.9749	2.9311	1.1702	-3.0145	-0.31147	0.4086	-3.0627	0.28336	3.741	0.33719	-0.05942	2.08E-15	
10	C163	5.3928	-2.3116	4.4967	0.32043	-0.98601	2.9897	1.3781	2.2271	-1.4572	-1.9265	2.1087	-0.49074	-4.14E-15	
11	C167	6.5384	-3.1193	-0.60883	2.4016	1.5316	-1.6575	0.49599	-0.92248	3.4164	-2.2906	-1.6113	0.67017	8.53E-16	
12	C144	-4.2809	-2.1855	-3.537	-0.05328	0.30189	1.0164	1.0039	-2.8895	-2.4345	-1.8862	0.94857	2.5747	7.85E-15	
13	C148	-2.6679	-2.3276	0.88478	-6.5421	0.67068	0.8113	0.51837	-0.91512	0.26234	-0.76997	-2.0472	-2.1479	3.30E-15	
14	C150	-3.774	-4.7669	-3.597	-0.45063	1.9031	1.1116	-1.5274	2.9032	1.9927	1.7823	1.9239	0.18468	8.20E-16	
15															

Note: The name of columns could be “PC1, PC2,...” or “Comp 1, Comp 2, ...”

Step 1: Prepare the Input Data

2. Prepare the metadata file with two columns (“**Sample**” and “**Group**”).

	A	B	C
1	Sample	Group	
2	C126	Group_1	
3	C132	Group_1	
4	C152	Group_1	
5	C153	Group_1	
6	C127	Group_2	
7	C135	Group_2	
8	C138	Group_2	
9	C160	Group_2	
10	C163	Group_2	
11	C167	Group_2	
12	C144	Group_3	
13	C148	Group_3	
14	C150	Group_3	
15			

Step 2: Upload the Data

1. Select tab “Scores Plot”.
2. Select subtab “Data Upload”.
3. Upload the metadata file and PCA/PLS-DA scores file.

Data Upload Plot Appearance

Axis Labels Axis Limits Axis Breaks

Save Plot

Upload Metadata File:

BROWSE... pca_metadata.csv

Upload complete

Upload PCA Score File:

BROWSE... pca_score.csv

Upload complete

Step 3: Modify the Scores Plot

Data Upload

Plot Appearance

Axis Labels

Axis Limits

Axis Breaks

Save Plot

Select and Order Group Levels

Group_1

Group_2

Group_3

Color for Group_1 Group:

#386CB0

Color for Group_2 Group:

#FDB462


Color for Group_3 Group:

#7FC97F

Legend label for Group_1 Group:

Group_1

Legend label for Group_2 Group:



Step 3: Modify the Scores Plot

1. Select subtab "Plot Appearance".
 - Select the option "Select and Order Group Levels" to modify the group order.

Data Upload

Plot Appearance

Axis Labels

Axis Limits

Axis Breaks

Save Plot

Select and Order Group Levels

Group_2

Group_1

Group_3

Color for Group_1 Group:

#386CB0

Color for Group_2 Group:

#FDB462


Color for Group_3 Group:

#7FC97F

Legend label for Group_1 Group:

Group_1

Legend label for Group_2 Group:



Step 3: Modify the Scores Plot

1. Select subtab "Plot Appearance".
 - Select the option "Color for Group" and select the appropriate color for each group.

[Data Upload](#)[Plot Appearance](#)

[Axis Labels](#)[Axis Limits](#)

[Axis Breaks](#)[Save Plot](#)

Select and Order Group Levels

Group_2Group_1Group_3

Color for Group_1 Group:

#369E19

Color for Group_2 Group:

#C73E1C


Color for Group_3 Group:

#787464

Legend label for Group_1 Group:

Group_1

Legend label for Group_2 Group:



Step 3: Modify the Scores Plot

1. Select subtab "Plot Appearance".
 - Select the option "Legend label for group" and modify the name for each group.

Data Upload

Plot Appearance

Axis Labels

Axis Limits

Axis Breaks

Save Plot

Select and Order Group Levels

Group_2

Group_1

Group_3

Color for Group_1 Group:

#369E19

Color for Group_2 Group:

#C73E1C

Color for Group_3 Group:

#787464

Legend label for Group_1 Group:


Non Survival

Legend label for Group_2 Group:

Survival

Legend label for Group_3 Group:

Control





Note: Legend label could contain the space.

Step 3: Modify the Scores Plot

1. Select subtab "Plot Appearance".

- Select or Remove select the option "Display 95% confidence ellipse" to display or remove the 95% confidence ellipse of groups.



Note: The 95% confidence ellipses are displayed as a default setting.

Step 3: Modify the Scores Plot

1. Select subtab "Plot Apperance".
 - Select option "Point Size" to resize the point size.
 - Select option "Plot Theme" to select the approriate theme.

Legend label for Group_3 Group:

Uninfected


☒ Display 95% confidence ellipse

Point Size:

5

Plot Theme:

theme_Publication



Step 3: Modify the Scores Plot

2. Select subtab “Axis Labels”.

Data Upload

Plot Appearance

Axis Labels

Axis Limits

Axis Breaks

Save Plot

X-axis Label:

The first component (A %)

Y-axis Label:

The second component (B %)

Axis Label Size:


15

☒ Axis bold

Tick Label Size:

15

☐ Tick bold



Step 3: Modify the Scores Plot

2. Select subtab “Axis Labels”.

- Input the % variance of PC1 and PC2 to “X-axis Label” and “Y-axis Label”, respectively.

[Data Upload](#)[Plot Appearance](#)

Axis Labels

[Axis Limits](#)[Axis Breaks](#)

[Save Plot](#)

X-axis Label:

The first component (26.4 %)

Y-axis Label:

The second component (15.2 %)

Axis Label Size:


15

☒ Axis bold

Tick Label Size:

15

☐ Tick bold



Note: The % variance of PC1 and PC2 were visualization output from e.g., MetaboAnalyst and NetworkAnalyst.

Step 3: Modify the Scores Plot

2. Select subtab “Axis Labels”.

- Modify the sizes of axis labels and tick labels in “Axis Label Size” and “Tick Label Size”, respectively.

[Data Upload](#)[Plot Appearance](#)

Axis Labels

[Axis Limits](#)[Axis Breaks](#)

[Save Plot](#)

X-axis Label:

The first component (26.4 %)

Y-axis Label:

The second component (15.2 %)

Axis Label Size:


20

☒ Axis bold

Tick Label Size:

20

☐ Tick bold



Step 3: Modify the Scores Plot

2. Select subtab “Axis Labels”.

- Modify the format (bold or not) of axis labels and tick labels at “Axis bold” and “Tick bold”, respectively.

Data Upload

Plot Appearance

Axis Labels

Axis Limits

Axis Breaks

Save Plot

X-axis Label:

The first component (26.4 %)

Y-axis Label:

The second component (15.2 %)

Axis Label Size:


20

☒ Axis bold

Tick Label Size:

20

☐ Tick bold



Note: The default setting of Axis label is bold and for the tick label size is not bold.

Step 3: Modify the Scores Plot

3. Select subtab “Axis Limit” to modify the limits of each axis.

[Data Upload](#)[Plot Appearance](#)

[Axis Labels](#)

Axis Limits

[Axis Breaks](#)


[Save Plot](#)

X-axis Minimum:

X-axis Maximum:

Y-axis Minimum:

Y-axis Maximum:



Step 3: Modify the Scores Plot

4. Select subtab “Axis Break” to modify the break of each column.

[Data Upload](#)[Plot Appearance](#)

[Axis Labels](#)[Axis Limits](#)[Axis Breaks](#)


[Save Plot](#)

X-axis Breaks:

6

Y-axis Breaks:

4



Step 4: Export the Scores Plot

1. Select subtab “Save Plot”.
2. Select the appropriate sizes (**Plot Width and Plot Height**) and resolution for the figure.

[Data Upload](#) [Plot Appearance](#)

[Axis Labels](#) [Axis Limits](#) [Axis Breaks](#)

Save Plot


Plot Width (in pixels):

Plot Height (in pixels):

DPI for Saving:

Format:

[↓ DOWNLOAD PLOT](#)



Note: Figure should be exported at 300 dpi or higher.

Step 4: Export the Scores Plot

3. Select figure format for the exported figure (.png, .tiff, .pdf, .svg).
4. Select option “Download Plot” to export the figure.

[Data Upload](#)[Plot Appearance](#)

[Axis Labels](#)[Axis Limits](#)[Axis Breaks](#)

Save Plot


Plot Width (in pixels):

Plot Height (in pixels):

DPI for Saving:

Format:

↓ DOWNLOAD PLOT



Note: The .svg format could be further processed with an appropriate toolkit.